

Restore to a file to a non-production library

Once you have located the backup needed for the restore create a temporary library

CRTLIB TEMP1 (you can call this anything you want)

You will restore to this library so you don't overwrite production data.

```
RSTOBJ OBJ(XXXX) SAVLIB(ABCXXXX) DEV(TAPXX) MBROPT(*ALL) ALWOBJDIF(*ALL)
RSTLIB(TEMP1)
```

OBJ=name of file you are restoring

SAVLIB=Original library name

DEV=the name of the device the backup is in

RSTLIB=the temporary library you created earlier

This parameter is very important so you don't overwrite your production data.